ENGLISH

LIGHT-CURE FLUORIDE RELEASING DENTAL RESTORATIVE MATERIAL

BEAUTIFIL II offers greatly improved handling properties and superb polishability. Other properties like fluoride release and recharge, wear resistance, outstanding radio-opacity, brilliant aesthetics because of controlled light diffusion, remain characteristics of BEAUTIFIL.

PRECAUTIONS FOR DENTAL PERSONNEL AND PATIENTS

1. Do not use this material on patients with known allergies to these types of chemicals. (This pertains to operators also.)
2. If any eruption or skin inflammation occurs on a patient when using this product, immediately discontinue use and have the patient seek medical advice. (The same precaution is applied for the operators)
3. Avoid contacting soft tissues, skin or eyes. In case of contact, rinse immediately with copious amounts of water and seek medical advice if necessary.

APPLICATION

Restoration of all anterior and posterior cavity classes I through V, wedge shaped defects and root caries.
Direct laminate veneer, Core build-up, Repair of restorations.

DIRECTIONS FOR USE

1. After cleaning tooth surface, take shade while the tooth is still moist, make cavity preparation, isolate prepared tooth with a rubber dam and protect the dentin close to the pulp with calcium hydroxide.
2. Apply bonding agent according to the manufacturer's instructions.
3. Syringe : Dispense the necessary amount of material from the syringe. After dispensing, immediately close the syringe cap then fill the cavity. Tips : Insert the selected shade tip into the Shofu TIPS dispensing gun or an other appropriate gun. Dispense slowly with consistent pressure.
4. If cavities exceeding 2.0mm depth, incremental placement no more than 2.0mm is recommended.
5. Light cure for 20 sec (halogen lamp) or 10 sec (high power type LED light). If the surface of the cavity is large, draw the light curing tip with a fanning action several times across the restoration.
6. Remove excess material with "Dura White Stones" or "T&F Hybrid Points". Final polishing is obtained by "CompoMaster".

NOTE IN USE

1. Use this product within the expiration date.
2. Keep away from direct sunlight and store at room temperature (1-30 °C).
3. If a cavity is deep, protect it with a pulp capping agent. (Do not use eugenol containing liner.)
4. Wear protective glasses when light curing.
5. Do not use this product in any other purpose except the purpose mentioned in the application.
6. This product is intended for dental operators use only.

COMPOSITION

Filler:	Multi-Functional glass filler and S-PRG filler based on fluoroboroaluminosilicate glass
Filler loading:	83.3wt% (68.6 vol%)
Particle size range:	0.01 to 4.0µm
Mean particle size:	0.8µm
Base resin:	Bis-GMA/TEGDMA resin

SHADES

A1, A2, A3, A3.5, A4, B1, B2, C2, Inc (Incisal), BW(Bleaching White), A20, A30

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A1, A2, A3, A3.5, A4, B1, B2, C2, Inc (incisal), BW (blanco blanqueador), A20, A30
A1, A2, A3, A3.5, A4, B1, B2, C2, Inc (Incisal), BW (Bleaching White), A20, A30