HC Primer

EN Instructions for use
ES Instrucciones de uso
RO Instrucțiuni de utilizare

DE Gebrauchsanleitung
IT Istruzioni per l'uso
SV Bruksanvisning

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Primer for CAD/CAM Hybrid Ceramic Restorations

PRODUCT DESCRIPTION

HC Primer is a single component primer developed to treat the bonding surface of CAD/CAM ceramic filler contained composite, to ensure a strong bond to abutments while luting with adhesive resin cement. It is recommended to be used in conjunction with resin cement system including self-etching primer for better results.

INDICATIONS FOR USE

Priming of CAD/CAM hybrid ceramic restorations

CONTRAINDICATIONS

Priming of CAD/CAM all-ceramic restorations

PRECAUTIONS AND WARNINGS

- Do not use this product on patients with known allergies to this material, methacrylate monomer and/or acetone.
- Operators with known allergies to this material and/or methacrylate monomer should exercise caution while using this product.
- In case of any inflammation or other allergic reactions in either patient or operator, immediately discontinue use and seek medical advice.
- Avoid inhaling vapor. If inhalation causes severe symptoms, including headache, dizziness, nausea or breathing difficulty, immediately seek medical attention.
- Use protective plastic gloves and glasses to avoid sensitization to this product. In case of accidental contact with oral soft tissue or skin, immediately blot with alcohol moistened cotton ball, and rinse with plenty of water. Medical (dental) gloves provided some protection against direct contact with this material; however some monomers may penetrate through the gloves very quickly to reach the skin. Should this material come in contact with the gloves, immediately remove them and wash hands with plenty of water. In case of accidental contact with eyes, immediately flush eyes with plenty of water and seek medical advice.
- Avoid any source of ignition since this product is flammable.
- Follow the instructions for use of each dental material, instrument or equipment used in conjunction with this product.
- Do not use this product for any purposes other than specifically outlined in these instructions for use.
- Use this product within the expiration date indicated on the package. (Example ♣ YYYY-MM-DD → Year-Month-Date of the expiration date)
- This product is intendeded for use by dental professionals only.

DIRECTIONS FOR USE

- Pretreatment of restoration
 - Sandblast the bonding surface of SHOFU BLOCK/DISK HC restorations with alumina sand (50-100 µm) at a pressure of 0.2-0.3 MPa. Clean with ultrasonic cleaning and air dry.
 - Or etch the bonding surface with 5 % hydrofluoric acid for 20 sec. Then, spray water, clean with ultrasonic cleaning and air dry.
 - Avoid using hydrofluoric acid in mouth.
 - Follow the respective CAD/CAM restorative manufacturer's instructions.
- Dispense 1-2 drops of the HC Primer into a dish. Apply a single, uniform coat of the Primer over the entire surface to be bonded with a brush. Then, gently air-dry the material surface until it no longer moives, while trying to spread it thinly and evenly to ensure optimal fitting of the restoration.
 - Allow the primer to reach room temperature before use when it is refrigerated, to avoid spewing of the Primer from the nozzle tip.
 - Tightly cap the bottle immediately after each use as it contains volatile agents.
 - Avoid contamination of the adherend surface of the restoration. In case of contamination with saliva or blood, rinse with water and thoroughly clean the surface with an alcohol moistened cotton ball. Then, dry and reapply the material.
- Light-cure using a dental light-curing unit. In the case where the adherend area is large, light-cure in segments.
 - Specifications of light-curing units that can be used for this material and their standard curing time are as follows.

Light-curing unit	Halogen	LED
Light-curing time	20 sec.	10 sec.
Wavelength	400-500 nm	440-490 nm
Light Intensity	≥ 500 mW/cm ²	≥ 1000 mW/cm ²

 - Operators must use protective glasses or light shielding plate to avoid direct exposure to the curing light.
- Placement of the restoration
 - Place the restoration onto the abutment or cavity following the instructions for use of the specific resin cement (eg. ResiCem). Then, light-cure each surface for the prescribed time.
 - Use a resin cement system with a primer for natural teeth.
 - Protect patients' eyes from the curing light with appropriate eye wear.

COMPOSITION

UDMA, MMA, Acetone, Polymerization initiator and others

STORAGE

Store in refrigerator (1-10 °C/34-50 °F). Avoid high humidity and keep away from direct sunlight and any source of ignition.

CAUTION: US Federal law restricts this device to sale by or on the order of a dental professional.

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