- EN** Instructions for use
DE Gebrauchsanleitung
FR Mode d'emploi
NL Gebruiksaanwijzing
ES Instrucciones de uso
IT Istruzioni per l'uso
RU Инструкция по применению

- PL** Instrukcja obsługi
RO Instrucțiuni de utilizare
SV Bruksanvisning
EL Οδηγίες χρήσης
PT Instruções de uso
ZH 使用说明书

05/2019-03
71803**EN Metal Bond Porcelain**

READ CAREFULLY BEFORE USE

INDICATIONS

High fusing porcelain for fabrication of metal bonded ceramic restorations

PRECAUTIONS

- If any inflammation or other allergic reactions occur on the operator while using this product, immediately discontinue use and seek medical advice.
- Avoid contact with skin or eyes. In case of accidental contact with eyes, immediately flush the eyes with plenty of water and seek medical advice.
- Use local dust extractor, dust protective mask, etc while grinding or polishing this product to avoid harmful influence of dust on the human body.
- This product is intended for use by dental professionals only.

NOTE ON USE

- Use this product in conjunction with the alloys having coefficient of thermal expansion of $13.6 \times 10^{-6} K^{-1}$ - $15.2 \times 10^{-6} K^{-1}$.
- The firing conditions may vary due to the different designs and operating voltages of the porcelain furnaces. It is essential to carry out test firing before using the porcelain for actual restoration.
- Second layer of Opaque must be applied evenly in a sufficient thickness to mask the color of the underlying metal coping.
- The Body build-up should be done to the final anatomical shape of the tooth in order to ensure the exact dimension for the Enamel porcelain.
- When large amount of porcelain is built up, set longer drying time at furnace entrance in order to avoid explosion of the porcelain from the inside due to remaining liquid components.
- Enamel porcelain should be built up slightly oversized in consideration of firing shrinkage.
- When the area requiring additional build-up is dry, wet the area with water using a brush.
- In the case where Opaque Porcelain is applied generally, set longer drying time at furnace entrance and preheating time in order to burn the liquid components thoroughly.
- Mix Opaque pastes sufficiently well with a plastic spatula before use.
- Use distilled water for Correction porcelains. Do not use VINTAGE MIXING LIQUID or MIXING LIQUID HC for Correction porcelains since they affect the shade of the restoration.
- Do not use VINTAGE MODELLING LIQUID since it affects the shade of the restoration.
- Wear protective glasses etc while grinding or polishing this product to avoid damage to eyes.
- Do not use this product for any purposes other than specifically outlined in the INDICATIONS in these instructions for use.

STORAGE

Store at room temperature (1-30 °C/34-86 °F). Avoid direct sunlight.

CAUTION: US Federal law restricts this device to sale by or on the order of a dental professional.

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


PHYSICAL PROPERTIES

Coefficient of thermal expansion Glass transition point

		Coefficient of thermal expansion (25–500°C)	Glass transition point (°C)
Base Opaque	2.Firing	11.5 × 10 ⁻⁶ K ⁻¹	598
	4.Firing		
Shade Opaque	2.Firing	12.4 × 10 ⁻⁶ K ⁻¹	594
	4.Firing		
Margin	2.Firing	13.0 × 10 ⁻⁶ K ⁻¹	594
	4.Firing		
Cervical, Opaque Dentin, Body, Enamel, Translucent, Uni-Layer	2.Firing	12.3 × 10 ⁻⁶ K ⁻¹	557
	4.Firing		
Cervical Trans (Cervical Trans, Gum)	2.Firing	12.9 × 10 ⁻⁶ K ⁻¹	548
	4.Firing		
Correction	2.Firing	11.7 × 10 ⁻⁶ K ⁻¹	528
	4.Firing		

FIRING SCHEDULE

	Drying and Preheating furnace entrance (min)	Temperature rising speed (°C/min)	Firing schedule
Base Opaque	5 to 7	50 to 60	500°C 940 to 960°C 2 min
Shade Opaque	5 to 7	50 to 60	500°C 920 to 940°C 1 to 2 min
Margin	5 to 7	50 to 60	650°C 920 to 940°C
Cervical	5 to 7	50 to 60	650°C 880 to 900°C
Body, Opaque Dentin, Opal, Opal Aqua, Opal Effect, Enamel Effect, Incisal, Uni-Layer, Color Effect, Cervical Trans	5 to 7	50 to 60	650°C 900 to 920°C
Self-glaze	5 to 7	50 to 60	650°C 880 to 900°C
Gum (First)	5 to 7	50 to 60	650°C 880 to 900°C
Gum-glaze	5 to 7	50 to 60	650°C 860 to 880°C
Correction (First)	5 to 7	50 to 60	450°C 750 to 770°C
Correction-glaze	5 to 7	50 to 60	450°C 750 to 770°C

SYSTEM AND SHADES

System	Shades
Base Opaque (1 shade)	BASE
Shade Opaque (21 shades)	A ₁ O, A ₂ O, A ₃ O, A _{3.5} O, A ₄ O, rootAO, B ₁ O, B ₂ O, B ₃ O, B ₄ O, C ₁ O, C ₂ O, C ₃ O, C ₄ O, D ₂ O, D ₃ O, D ₄ O, W ₀ O, W ₁ O, W ₂ O, W ₃ O
Opaque Effect (7 shades)	OM-G, OM-O, OM-Br, OM-Y, OM-BG, OM-LP, OM-DP
Cervical (4 shades)	AC, BC, CC, DC
Margin (14 shades)	W ₀ M, NM, A ₁ M, A ₂ M, A ₃ M, A _{3.5} M, A ₄ M, rootAM, B ₂ M, B ₄ M, C ₂ M, C ₄ M, D ₃ M, CLM
Margin Effect (6 shades)	LPM, DPM, OrM, BrM, MLM, MDM
Body (21 shades)	A ₁ B, A ₂ B, A ₃ B, A _{3.5} B, A ₄ B, rootAB, B ₁ B, B ₂ B, B ₃ B, B ₄ B, C ₁ B, C ₂ B, C ₃ B, C ₄ B, D ₂ B, D