



SHOFU DENTAL CORPORATION  
1225 Stone Drive, San Marcos, California 92078-4059 U.S.A.  
Tel: (760) 736-3277 • Fax: (760) 736-3276

## **MonoCem™**

1. Compressive Strength: 44,500 p.s.i. / 310 MPa.
2. Film thickness: translucent shade 11 microns  
bleach shade 12.5 microns
3. Filler: 60%
4. Retention value (non-threaded titanium post): 32 kg
5. Flouride Release 7.9 ppm
6. Recommended for Zirconia and other reinforced all ceramics, Cerec, gold, in addition to pfm's
7. Number of total steps in procedure 3 (versus 7 for Unicem and 9 for Panavia)  
(Dentistry Today, March 2005)
8. Minimal waste (Ave. 0.2 ml/tip – CRA)
9. 100% polymerization