



JONATHAN EAGLE, DDS

Dr. Eagle is an active clinician in Grand Rapids, MI. On graduating from the University of Michigan School of Dentistry in 1990, he received training in Advanced Anterior Aesthetics and as a Master Aesthetic Dentist from the Las Vegas Institute. Dr. Eagle is certified to perform Invisalign treatment and received extensive training under Dr. Richard Litt, former director of the Orthodontics Department at the University of California at Berkeley and the University of Detroit. He also obtained advanced hands-on endodontic training from Dr. Stephen Buchanan at the University of Detroit. Dr. Eagle is a member of the West Michigan Dental Association, the Michigan Dental Association, the American Dental Association, and the American Academy of Cosmetic Dentistry.



FOR FREE INFORMATION:

800.921.4806
ext. 136
www.dps.
li/a/58N-405
Circle 405
on the card

EyeSpecial C-II

SHOFU DENTAL

A smart, digital dental camera, EyeSpecial C-II is designed to achieve predictable and consistent clinical photographs for case documentation, diagnosis and treatment planning, patient communication, and dental lab collaboration. Equipped with a 12 mega-pixel sensor and cutting-edge FlashMatic module—a proprietary system of ring and dual-point flashes—this camera demonstrates true-color reproduction with an excellent depth-of-field range. With its ultra-lightweight body design and large LCD touch-screen, EyeSpecial C-II allows clinicians to easily transition to 8 preset dental modes with one touch. Dr. Jonathan Eagle adopted the technology earlier this year and shares his experience with us.

DPS **When did you start using the EyeSpecial C-II in your practice?**

I purchased it in January. My motivation for getting it was so that I had a camera that was easier and simpler to use so my hygienist could document photographs for patients and get prior authorizations for treatments, rather than having me come down and do it for her. We were looking for a simpler and easier-to-use camera.

DPS **Had you used digital dental cameras in the past?**

Yes, I had used a digital SLR for years and years. We're very big about documenting cases as we go. We document before treating the tooth, we document when we've removed the previous restoration, and we document decay. It speeds the pathway for billing purposes with insurance and provides good documentation if there



was ever an issue. I also have an emphasis in cosmetic dentistry, so we obviously like to have some before-and-after photographs for those cases as well.

DPS **What features attracted you to EyeSpecial C-II?**

The biggest selling point for me is that it's much easier to use one-handed than a digital SLR.

You have the weight of the SLR, and you have to put your eye up to the viewfinder on SLR. When I was using an SLR, I was using it with a manual focus set in one position, and I had to move the camera in and out to gain focus. The auto-focus on EyeSpecial C-II is extremely accurate and intuitive to use, and the body of the camera is ultra lightweight. When you're trying to maneuver