

# Beautiful® II

## PINK & WHITE

FLUORIDE RELEASING, BIOACTIVE,  
NANO-HYBRID COMPOSITE



pure aesthetics



- GIOMER Technology—clinically proven to release and recharge fluoride, reduce plaque and neutralize acids
- Easy to place, sculpt and polish
- High flexural and compressive strength with superior wear resistance
- Natural fluorescence and high radiopacity, 3.4 al:mm

gingiva & enamel  
shade composite

Visit [www.shofu.com](http://www.shofu.com) or call 800.827.4638

# PINK & WHITE



Photo courtesy of Frank J. Milnar, DDS, AAACD

- Aesthetic gingival reproduction—  
Recession, exposed abutments and PFM crown margins, root erosion and lost papilla
- Shades are designed to meet the needs of all ethnicities
- Aesthetic solutions without surgery
- Perfect for Class V restorations

	G-Dark Pink	G-Light Pink	G-Brown	G-Orange	G-Violet
2.5g					
PN	Y2257	Y2256	Y2259	Y2260	Y2258

## PINK Prep Kit\*




\*Developed in cooperation with Frank J. Milnar, DDS, AAACD

- Designed to effectively prepare the surface for functional aesthetic results
- Five burs to create the ideal Minimally Invasive prep design



Photo courtesy of Frank J. Milnar, DDS, AAACD

- Achieve highly aesthetic direct restorations with minimal shade selections
- Create natural aesthetics with a simplified method for shade modification
- Modeled after the porcelain color concept—easy color matching
- Direct veneers made easy



	High-Value Translucent (HVT)	Translucent (T)	Low-Value Translucent (LVT)	Amber (AM)
2.5g				
PN	Y2251	Y2250	Y2252	Y2253

# Beautiful® II PINK & WHITE

### Beautiful II PINK Kit

Part Number **GMK1**

- 5 (2.5g) Beautiful II Gingiva syringes
- 1 (2g) Beautiful Flow Gingiva F02 syringe

#### Bonus Items Included!

- 4 Beautiful Flow Plus F00 & F03 A2 Tips
- 2 Beautiful II A2 Tips

### COMPLEMENTARY PRODUCT

#### PINK Prep Kit

Part Number **8205K1**

- 2 Minimally Invasive FG Robot Diamonds
- 3 FG Robot Diamonds
- 5-Hole Ruby Pink BurButler


### Beautiful II PINK & WHITE Kit

Part Number **GMK3**

- 5 (2.5g) Beautiful II Gingiva syringes
- 1 (2g) Beautiful Flow Gingiva F02 syringe
- 4 (2.5g) Beautiful II Enamel syringes

#### Bonus Items Included!

- 10 Beautiful Flow Plus F00 & F03 A2 Tips
- 10 BeautiBond Unite Dose
- 10 Super-Snap Singles Mini



Shofu Dental Corporation  
1225 Stone Drive  
San Marcos, CA  
92078-4059, USA