

What's Hot and What's Getting Hotter

BY HOWARD S. GLAZER, DDS, FAGD

MATERIAL

Beautiful® Flow Plus
Shofu Dental Corporation
1225 Stone Drive
San Marcos, CA 92078
800.827.4638
www.shofu.com

I have often said with regard to flowable resins that I wanted a flowable that was not a "runnable." Shofu has created just such a product with its new, non-flow, sculptable, Beautiful® Flow Plus composite resin. This is a great material as a final restoration for minimally invasive procedures in both the posterior and anterior dentition. You do not need to place another layer of composite resin on top. I am quite pleased with how well this material handles, especially when restoring root caries and the more conventional Class V lesions (whether caused by decay or abrasion). The material stays where it is placed and is available in two formulations: a non-flow sculptable (F00) or a low-flow (F03). There are nine shades for F00, including opaque, incisal, and bleach white, and 12 shades for F03, including a unique A0.5, milky, and cervical shades. I am quite impressed by the great diverse quality of these resins. This material also is useful in creating

cuspal height and marginal ridges. Beautiful Flow Plus contains the proprietary GIOMER chemistry and S-PRG fillers, which release and recharge fluoride like a glass ionomer. Both formulations will polish to a very high gloss, which is expected from the people who produce the highly acclaimed One Gloss and Super Snap Single polishing systems. This is a very useful resin product that allows us to truly emphasize the artistic element of the "art and science" of dentistry.



Reprinted with permission from AGD Impact, November 2010. On the web at www.agd.org.
© Academy of General Dentistry. All rights reserved. Foster Printing Service: 866-879-9144, www.marketingreprints.com.



SHOFU DENTAL CORPORATION
Call for Information 800 827 4638 www.shofu.com
SAN MARCOS CALIFORNIA


FINALLY, AN INJECTABLE HYBRID RESTORATIVE FOR ALL INDICATIONS

BEAUTIFIL Flow Plus

- Indicated in All Classes
- Physical Properties Rival Hybrids
- Precise Handling Stays Put
- Self-polishing "Leveling-effect"

BEAUTIFIL Flow Plus is the next step in the evolution of restorative materials, combining the delivery of a flowable and the strength, durability, and aesthetics equal to or better than leading hybrid composites.

Visit www.shofu.com or call 800.827.4638



Shofu Dental Corporation • San Marcos, CA

**TRIAL OFFER \$99.95
A \$160.00 VALUE!**

