


# Glaslonomer FX-II Capsule

Enhanced Direct Restorative



# Glaslonomer FX-II Capsule

Enhanced Direct Restorative

An aesthetic all-purpose glass ionomer cement for direct restoration, designed for Minimal Invasive dentistry.

- **No conditioner** required with equal bond strength to dentin and enamel
- Thermal expansion similar to natural teeth
- Ideal handling and packability
- Ready to use immediately after mixing – no wait
- Radiopaque and low solubility
- Increased fluoride release and recharge by 39% over the market leader
- Highest value in hardness and equal compressive strength when compared to leading competitors
- 3 year shelf life
- Significantly decreases acid erosion levels

The Glaslonomer FX-II Capsule is an all-purpose enhanced direct restorative indicated for non-load bearing Class I & Class II restorations in permanent dentition and Class III & Class V cavities with wedge shaped defects. Ideal for pediatric and geriatric patients as a core build-up and base/liner material.


Visit [www.shofu.com](http://www.shofu.com) or call **800.827.4638**



One: Before



Two: After



Fluorescence

Photos courtesy of George Freedman DDS and Shofu, Inc.

CASE STUDY



#### Order Information

50 Pre-mixed Capsules

Shades	Item No.
A2	3270
A3	3271
A3.5	3272
B2	3273



Shofu Dental Corporation  
1225 Stone Drive  
San Marcos, CA  
92078-4059, USA

SNF881-1213

Enhanced Direct Restorative