


OneGloss®

ONE SYSTEM—ONE STEP
FINISHING AND POLISHING



EASY PROCEDURE



Finish and polish in one step
Does not scratch enamel



Removes cement or composite flash
Removes intrinsic stain



Does not crumble
Disposable or reusable shank

OneGloss®

OneGloss is an aluminum oxide one-step finisher and polisher that provides the ideal finish for all types of composite and cemented restorations simply by altering the contact pressure as you go from finishing to polishing. Available in three different shapes and durable enough to withstand polishing an entire arch.

OneGloss

- Durable and disposable
- Complete kit with stainless steel CA mandrels (no metal hub)




OneGloss PS

- Durable, disposable and individually wrapped
- Mounted on sturdy plastic mandrels (no metal hub)



Visit www.shofu.com or call **800.827.4638**

SHAPE CATEGORY		Cup	Midi-Point	IC	Mini-Point	CA Mandrel
Shape						
ISO #		060	060	100	030	
OneGloss						
CA (50 pk)	PN	0181	0182	0183	0186 (12 pk)	0435 (6 pk)
OneGloss PS						
CA (50 pk)	PN	0176	0177	0178	-	-



Order Information

- Part Number 0180 – OneGloss Set
25 Cups, 25 Midi-Points, 25 IC and 3 CA mandrels
- Part Number 0175 – OneGloss PS Assortment Kit
15 Cups, 15 IC and 20 Midi-Points
- Part Number 0185 – OneGloss SuperBuff Combo Set
20 OneGloss Cups, Midi-Points and IC, 25 SuperBuff disks and 5 CA mandrels
- See chart above for individual shapes



Shofu Dental Corporation
1225 Stone Drive
San Marcos, CA
92078-4059, USA

SN69-0610

Finishing and Polishing