



ABRASIVE QUICK GUIDE - CHAIRSIDE

OPERATORY	CONTOURING/FINISHING	POLISHING	HIGH-LUSTER POLISHING
ENAMEL			
	PINK Prep Kit	OneGloss Dura-White CeraMaster Coarse Ceramisté Standard, Ultra	DirectDia Polishing Paste with Super-Snap Buff Disk Super-Snap SuperBuff Disk CeraMaster Ceramisté Ultra II
FILLINGS			
COMPOSITES	Dura-Green/Dura-Green Fine Dura-White OneGloss Super-Snap black/violet Robot Carbide 8/16 blades T&F Hybrid Points CompoSite Artistic Composite Kit	OneGloss Super-Snap green Super-Snap X-Treme green Robot Carbide CTF 30 blades CompoMaster Coarse CompoSite Fine	Super-Snap red Super-Snap X-Treme red CompoSite Polishing Paste with Super-Snap Buff Disk CompoMaster SuperBuff
AMALGAMS	Dura-Green	Brownie/Greenie	Supergreenie
GLASS-IONOMERS	Dura-White T&F Hybrid Points Robot Carbides		Ceramisté Ultra II
INLAYS, CROWNS, BRIDGES, FRAMEWORKS			
TITANIUM	Dura-Green Hardie Black	Brownie/Greenie Hardie Brown	Supergreenie
PRECIOUS ALLOYS	Dura-Green	Brownie/Greenie	Supergreenie
NON-PRECIOUS ALLOYS	Dura-Green Hardie Black	Brownie/Greenie Hardie Brown	Supergreenie
VENEERS, FACINGS, INLAYS, CROWNS			
CERAMICS (Zirconia, Lithium Disilicate, **IPS e.max®, **Bruxzir®, Alumina)	Dura-Green DIA Dura-Green Dura-White ZiLMaster	CeraMaster Coarse ZiLMaster	CeraMaster Ceramisté Ultra II DirectDia Polishing Paste ZiLMaster
PORCELAINS	Dura-Green DIA Dura-Green/Dura-Green Fine Dura-White T&F Hybrid Points	CeraMaster Coarse Ceramisté Standard, Ultra Soft PA & PB	CeraMaster Ceramisté Ultra II DirectDia Polishing Paste
COMPOSITES & HYBRID CERAMICS	Dura-Green/Dura-White T&F Hybrid Points OneGloss	CompoMaster Coarse CompoSite/CompoSite Fine Ceramisté Soft PA & PB OneGloss	CompoMaster CompoSite Polishing Paste with Super-Snap Buff Disk CompoSite Fine
LAB COMPOSITES (*Lava™ Ultimate, Ceramage)	Dura-Green DIA Dura-Green	Ceramisté Soft PB CompoMaster Coarse *Dura-Polish	CompoMaster *Dura-Polish DIA
PARTIAL/FULL DENTURES			
DENTURE ACRYLICS	*Robot Carbide HP Cutters *AcryPoint, dark grey and brown	*AcryPoint, light grey	
MOUTHGUARDS/SPLINTS			
ACRYLICS (*Under Armour®)	Dura-Green DIA Dura-Green *AcryPoint, dark grey	*AcryPoint, brown *AcryPoint, light grey	

■ one-step procedure

■ one-step procedure and disposable

* Not for intraoral

**All trademarks belong to their manufacturers

■ two-step procedure

■ three-step procedure

■ multi-step procedure