

ABRASIVES

ABRASIVE QUICK GUIDE - LABORATORY

LABORATORY	CONTOURING/FINISHING	POLISHING	HIGH-LUSTER POLISHING
VENEERS, FACINGS, INLAYS, CROWNS			
CERAMICS (Zirconia, Lithium Disilicate, **IPS e.max®, **Bruxzir®, Alumina)	Dura-Green DIA ZiLMaster	Dura-Polish ZiLMaster	Dura-Polish DIA ZiLMaster
PORCELAINS	Dura-Green DIA Dura-Green/Dura-Green Fine Dura-White Robot Carbide CTF T&F Hybrid Points	CeraMaster Coarse CeramiSté Standard, Ultra Soft PA & PB Dura-Polish	CeraMaster CeramiSté Ultra II DirectDia Polishing Paste Dura-Polish DIA
LAB COMPOSITES (**Lava™ Ultimate, Ceramage) HYBRID CERAMICS (Shofu HC)	Dura-Green Dura-White OneGloss Robot Carbide CTF T&F Hybrid Points Robot Carbide HP Cutters	CompoMaster Coarse CeramiSté Soft PA & PB OneGloss Dura-Polish	CompoMaster CompoSite Paste Dura-Polish DIA
INLAYS, CROWNS, BRIDGES, FRAMEWORKS			
PRECIOUS ALLOYS	Lab Pink Lab White	Brownie/Greenie	Supergreenie Dura-Polish
NON-PRECIOUS ALLOYS	Robot Carbide Lab Coral Hardie Black	Hardie Brown	Dura-Polish
PARTIAL DENTURES			
MODEL CASTING	Robot Carbide Lab Coral Hardie Black	Hardie Brown	
PARTIAL/FULL DENTURES			
DENTURE ACRYLICS	Robot Carbide HP Cutters AcryPoint, dark grey and brown	AcryPoint, light grey	Dura-Polish
ORTHO APPLIANCES/SPLINTS/MOUTHGUARDS			
ACRYLICS (**Under Armour®)	AcryPoint, dark grey Robot Carbide HP Cutters	AcryPoint, brown AcryPoint, light grey	

■ one-step procedure

■ one-step procedure and disposable

**All trademarks belong to their manufacturers

■ two-step procedure

■ three-step procedure