



ABRASIVES

ACRYPOINT LAB POLISHERS

SHAPE CATEGORY	BP1	PC2
SHAPE		
ISO#	100	050
DARK GREY (coarse)		
HP (6 PK)	PN 0426	-
HP (12 PK)	PN -	0423
BROWN (medium)		
HP (6 PK)	PN 0427	-
HP (12 PK)	PN -	0424
LIGHT GREY (fine)		
HP (6 PK)	PN 0428	-
HP (12 PK)	PN -	0425

- Available in 3 silicon carbide grits
- Use **Dark Grey (coarse)** for adjustments, **Brown (medium)** for finishing & **Light Grey (fine)** for polishing
- Use **AcryPoint for:** acrylic resin denture material

ACRYPOINT LAB POLISHERS

PN

Product

See chartHP Shank (6 per pack)

BP1

See chartHP Shank (12 per pack)

PC2