

# Beautiful-Bulk Flowable

## Bulk Fill Flowable Composite with Bioactive GIOMER Technology



### Beautiful-Bulk Flowable

Bulk Fill Flowable Composite with Bioactive GIOMER Technology


Shofu is dedicated to manufacturing the highest quality products utilizing our innovative technologies to enhance minimally invasive cosmetic dentistry.

shofu.com  
800.827.4638  
customer-service@shofu.com




### FUNCTIONAL

**Beautiful-Bulk Flowable** is the ideal light-cure dentin replacement material for deep restorations. For years glass ionomer was used for dentin replacement, but long curing times caused significantly reduced productivity, making this impractical until resin-based alternatives emerged. With physical properties similar to dentin, 4mm depth of cure, and the lowest shrinkage stress in this category, **Beautiful-Bulk Flowable** is the perfect choice for sensitivity-free dentin replacement and maximum clinical efficiency.




### BIOACTIVE GIOMER TECHNOLOGY INSIDE

Shofu's proprietary Giomer Technology gives **Beautiful-Bulk Flowable** anti-bacterial and acid neutralizing effects to minimize problems caused by bacterial proliferation such as secondary decay and restorative failure. A Univ. of Florida Dental School clinical study showed no secondary decay or failures after 8 years! **Beautiful-Bulk Flowable** is the only bulk flow product to combine bioactive benefits with speed and clinical efficiency.




Minimal Bacteria (Asahi University)




Advanced Bacterial Colonization

#### Fluoride Release & Recharge from GIOMER



\*1. Gordan V.V., et al. A clinical evaluation of a self-etching primer and a Giomer restorative material: Results at 8 years. JAM Dent. Assoc. 2007;138:621-627

## Beautiful-Bulk Flowable



### Beautiful-Bulk Flowable

2.4 g syringe

Shade

Universal

Dentin

Part Number

2030

2031

0.23 g x 20 tips

Shade

Universal

Dentin

Part Number

2028

2029



Scan Here  
Instructions for Use

# Beautiful-Bulk Restorative

## Bulk Fill Posterior Composite with Bioactive GIOMER Technology



### Beautiful-Bulk Restorative

Bulk Fill Posterior Composite with Bioactive GIOMER Technology

Shofu is dedicated to manufacturing the highest quality products utilizing our innovative technologies to enhance minimally invasive cosmetic dentistry.

shofu.com  
800.827.4638  
customer-service@shofu.com



### FUNCTIONAL


**Beautiful-Bulk Restorative** is the ideal metal-free material for bulk fill posterior restorations. This highly filled material (87.0%) offers low shrinkage (<1.7%) and wear resistance, 4mm depth of cure, as well as the combination of superior compressive and flexural strength for maximum durability. In addition, **Beautiful-Bulk Restorative's** esthetics are unmatched by any competitor.




Jack Griffin, DDS

### BIOACTIVE GIOMER TECHNOLOGY INSIDE

Shofu's proprietary GIOMER Technology gives **Beautiful-Bulk Restorative** long-term ion exchange for anti-bacterial and acid neutralizing effects to minimize problems caused by plaque and biofilm which may lead to restorative failure. A Univ. of Florida Dental School clinical study showed no secondary decay or failures after 8 years! **Beautiful-Bulk Restorative** is the only product in this market segment that gives the bioactive and therapeutic benefits of GIOMER Technology and the clinically efficiency of posterior bulk filling.




Minimal Bacteria (Asahi University)




Advanced Bacterial Colonization

#### Fluoride Release & Recharge from GIOMER



\*1. Gordan V.V., et al. A clinical evaluation of a self-etching primer and a Giomer restorative material: Results at 8 years. JAM Dent. Assoc. 2007;138:621-627

## Beautiful-Bulk Restorative



### Beautiful-Bulk Restorative

4.5 g syringe

Shade  
Universal  
A

Part Number  
 2034  
 2035

0.25 g x 20 tips

Shade  
Universal  
A

Part Number  
 2032  
 2033



Scan Here  
Instructions for Use