



SHOFU ROBOT® CARBIDE CTF-FG FINISHERS



- High quality tungsten carbide finishing head on a stainless steel shank
- Trim & finish subgingival cosmetic restorations without damaging the subgingival tissue
- Available in three blade types: 8, 16 (yellow-banded shank) & 30 (white-banded shank)
- 8 blades: rapid trimming & contouring; 16 blades: Fine Finishing (F); 30 blades: Ultra Fine Finishing (UF)
- 30 blades: Ultra Fine Polishing (UF)
- **Use Robot Carbide CTF-FG Finishers for:** cosmetic restorations

SHAPE CATEGORY	SH 132 F	SH 132 UF	SH 133	SH 133 UF
SHAPE				
HEAD SIZE (MM)	008	008	010	010
HEAD LENGTH (MM)	3.0	3.0	4.2	4.2
# OF BLADES	16	30	8	30
MAXIMUM RPM	300,000	300,000	450,000	450,000
FG (3 PK)	PN 0762-3	0763-3	0764-3	0766-3

SHAPE CATEGORY	SH 134	SH 134 F	SH 134 UF	SH 135	SH 135 F	SH 135 UF
SHAPE						
HEAD SIZE (MM)	014	014	014	014	014	014
HEAD LENGTH (MM)	6.0	6.0	6.0	9.0	9.0	9.0
# OF BLADES	8	16	30	8	16	30
MAXIMUM RPM	450,000	450,000	450,000	300,000	300,000	300,000
FG (3 PK)	PN 0767-3	0768-3	0769-3	0770-3	0771-3	0772-3

PRICING - SHOFU ROBOT® CARBIDE CTF-FG FINISHERS

PN

Product

See chartsCTF-FG Regular - 8 blades (3 per pack)

See chartsCTF-FG Fine - 16 blades (3 per pack)

See chartsCTF-FG Ultrafine - 30 blades (3 per pack)