


ABRASIVES

CERAMISTÉ®/CERAMISTÉ® SOFT POLISHERS



- Silicon carbide impregnated polishers
- Three-step system
- Use Ceramisté Standard for pre-polishing, Ultra (yellow band) for polishing & Ultra II (white band) for super-polishing
- Use Ceramisté Soft PA (green band) & Ceramisté Soft PB (blue band) for smoothing rough surfaces & fragile porcelain margins
- Use **Ceramisté for:** porcelains & enamel

SHAPE CATEGORY	Bullet	Cup	Knife	Wheel	Midi-Pt.	KN7	PC2	WH6	Unm. Wheel
SHAPE									
ISO#	040	065	090	090	060	125	050	120	220
STANDARD (no band)									
HP (12 PK)	PN	-	-	-	-	0261	0262	0263	0451
CA (12 PK)	PN	0251	251A	0252	0253	253A	-	-	-
FG (12 PK)	PN	-	251B	-	-	253B	-	-	-
ULTRA (yellow band)									
HP (12 PK)	PN	-	-	-	-	0264	0265	0266	0452
CA (12 PK)	PN	0254	254A	0255	0256	256A	-	-	-
FG (12 PK)	PN	-	254B	-	-	256B	-	-	-
ULTRA II (white band)									
HP (12 PK)	PN	-	-	-	-	0267	0268	0269	452A
CA (12 PK)	PN	0257	257A	0258	0259	259A	-	-	-
FG (12 PK)	PN	-	257B	-	-	259B	-	-	-
SOFT PA (green band)									
HP (12 PK)	PN	-	-	-	-	-	-	0564	0595
SOFT PB (blue band)									
HP (12 PK)	PN	-	-	-	-	-	0563	0565	0596

CERAMISTÉ®/CERAMISTÉ® SOFT POLISHERS

PN

Product

See chart. HP, CA, FG Shank (12 per pack)

See chart. Unmounted Wheel (12 per pack)

See chart. Soft PA/PB (12 per pack)

See chart. Soft PA/PB Unmounted Wheel (12 per pack)