









ABRASIVES

COMPOSITE/COMPOSITE FINE POLISHERS



- Use on microfilled & hybrid composites
- Remove orthodontic adhesives
- Two-step system to pre-polish & polish
- Use **CompoSite for:** composites

SHAPE CATEGORY		Bullet	Knife	Wheel	Disk	Cup	Pear-Point	IC	Midi-Point
SHAPE									
ISO#		040	090	090	160	080	033	-	040
REGULAR (no band)									
CA (12 PK)	PN	0191	0192	0193	0194	0195	0196	0197	0198
FG (12 PK)	PN	-	-	-	-	-	196B	-	0199
FINE (white band)									
CA (12 PK)	PN	0291	0292	0293	0294	0295	0296	0297	0298
FG (12 PK)	PN	-	-	-	-	-	296B	-	0299

COMPOSITE/COMPOSITE FINE POLISHERS

PN Product

See chart CA & FG Shank (12 per pack)