

## CORESHADE GLASIONOMER (CORE BUILD-UP BASE CEMENT)



- Supportive and protective base under stress-bearing fillings & restorations
- High compressive strength (147 MPa)
- Excellent bonding to tooth structure, metal pins & posts
- Gray in color for easy identification during placement & finishing
- High fluoride release & high level of radiopacity
- Biocompatible

### CORESHADE GLASIONOMER

#### PN

#### Product

**1115 . . . . .CoreShade Base Cement Kit**

25g gray powder, 10ml liquid, mixing pad, dispenser & instructions

**1116 . . . . .Gray Powder 25g**

**1117 . . . . .CoreShade Liquid 10ml**