

DentaLaze

High-Quality, Low Cost Wireless Soft-Tissue Diode Laser



DentaLaze

High-Quality, Low Cost
Wireless Soft-Tissue
Diode Laser

**SERVE
PATIENTS
INCREASE
PRODUCTION**

**GO HOME
ON TIME!**

Shofu is dedicated to manufacturing the highest quality products utilizing our innovative technologies to enhance minimally invasive cosmetic dentistry.

shofu.com
800.827.4638
customer-service@shofu.com

CLINICAL EFFICIENCY

Perform a variety of soft-tissue procedures with greater patient comfort and less healing time. The economically priced DentaLaze streamlines your workflow by being indicated for a multitude of Dental applications - both Clinical and Hygiene.* Featuring an output of up to 3W continuous power or 5W pulsed power at 808nm, this laser is sure to be a workhorse around your practice. The extended battery life enables the DentaLaze to run for up to 45 minutes on a single charge.

*Consult your state laws about regulations surrounding the use of lasers in your practice.



EASE OF USE

The wireless, aluminum DentaLaze is sturdy yet light weight and is controlled through an intuitive push button interface. The ergonomic design makes it perfect for either anterior or posterior applications. Its distinctive green aiming beam guarantees better visibility than competitors who feature a red aiming beam, which can be harder to see against a background of intraoral tissue or blood. The handheld, compact design allows the DentaLaze to be transported from room-to-room without the need for cords or foot pedals - and the single patient use fiber tips reduce the potential for cross-contamination.



ECONOMICAL PRICE

Laser Dentistry doesn't have to be expensive. DentaLaze costs significantly less than many other soft tissue diode laser competitors, yet has features which makes it more versatile. Surveys show that patients have a positive perception of lasers in the Dental practice.* This can lead to increased case acceptance, broader scope of treatment, shorter procedure time, and faster recovery time.

*Wigdor H. Patients' perception of lasers in dentistry. Lasers Surg Med. 1997;20(1):47-50. doi: 10.1002/(sici)1096-9101(1997)20:1<47::aid-lsm7>3.0.co;2-z. PMID: 9041507.

POPULAR PROCEDURES

For laser procedures, including:

- Hemostasis
- Canker Sores/Aphthous Ulcers
- Frenectomy
- Periodontal Sulcular Debridement
- Gingivectomy
- Gingival Troughing for Impression
- Incision/Excision
- Desensitization
- Root Canal Therapy
- And many other soft tissue procedures!

Online technique
videos are
available at
shofu.com/dentalaze

TECHNICAL SPECIFICATIONS

WAVELENGTH (±10NM)	DENTALAZE	808NM
MAX. POWER (±20%; WATT)	DENTALAZE	CW 3 WATT
		PULSE 5 WATT (PULSE DURATION 0.225 SECONDS; PULSE FREQUENCY 2.2HZ)
DUTY CYCLE		CW 100%; PULSED 50%
BATTERY PERFORMANCE		OUTPUT: 3.7V NOMINAL AND 1.4AH RUNNING TIME ~ 45 MINUTES (1.5W) 90 MINS. FULL RECHARGING TIME EXP. LIFE: 500 FULL CHARGE/DISCHARGE CYCLES
DIMENSION (MM, WITHOUT FIBER TIP)		LENGTH: 163; WIDTH: 22; HEIGHT: 22
CHARGING DOCK DIMENSIONS (MM)		HEIGHT: 45; DIAMETER:44
WEIGHT		HANDPIECE: 100 GM
LASER CLASSIFICATION		4
LIMITED WARRANTY		TWO-YEAR WARRANTY AGAINST MANUFACTURER DEFECTS
		ONE-YEAR BATTERY AND CHARGING STAND



shofu.com
800.827.4638
customer-service@shofu.com



DentaLaze

The DentaLaze Diode Laser Kit is available for only \$3,495

Contact kayceh@shofu.com
or call 800-827-4638
for more information



Scan Here
Instructions for Use



DENTALAZE DIODE LASER KIT (\$3,495)

Item Number: DL001

Included with the kit:

- 1 Laser Handpiece - (Main Control & Battery)
- 2 Goggles
- 1 Laser Shield
- 5 Fiber Optic Tips
- 50 Barrier Sleeves
- 1 Charging Stand
- 1 Power Adapter

FIBER OPTIC TIPS 25/PK (\$159)

Item Number: DL002

BARRIER SLEEVES 500/PK (\$99)

Item Number: DL003

GOGGLES 1/PK (\$249)

Item Number: DL004

LASER SHIELD (\$199)

Item Number: DL005

LASER HANDPIECE MAIN CONTROL BATTERY (\$249)

Item Number: DL006

CHARGING STAND (\$199)

Item Number: DL007

POWER ADAPTER (\$99)

Item Number: DL008

ACADEMY OF LASER DENTISTRY LASER SAFETY COURSE (\$150)

*Plus tax and shipping charges where applicable

TO ORDER A DENTALAZE FILL OUT THE FORM AND EMAIL TO KAYCEH@SHOFU.COM

Name _____ Date _____

Address _____ State _____ Zip _____

Phone _____ Alt. Phone _____ Email _____

Amex/Visa/MC Account # _____ Exp. Date _____ Val. Code _____

Subtotal (Tax & Shipping not included) _____

Customer Signature _____ Date _____