

ABRASIVES

DURA-GREEN® DIA STONES



- NEW FG SHAPES!
- Vitrified bonded diamond abrasive
- Improved cutting, better than carborundums
- Reduced chipping & contamination
- Use **Dura-Green DIA Stones** for: contouring zirconia, porcelain & ceramics

SHAPE CATEGORY	CY5	FL3	IC5	KN7	RD2	RD3	TC2	TC4	WH6
SHAPE									
ISO#	056	030	070	125	024	036	025	025	130
HP (2 PK*/3 PK)	PN 0156*	0158	0162	0154*	-	0159	0161	0155	0157*
CA (3 PK)	PN -	0163	-	-	-	0164	0165	-	-
FG (3 PK)	PN -	0168	-	-	0167	-	-	-	-

DURA-GREEN® DIA STONES

PN

Product

See chartsHP, CA, FG Shank

0153HP Shank Dura-Green DIA Assortment
CY5, KN7, TC4

0166CA Shank Dura-Green DIA Assortment
FL3, RD3, TC2

0504Shofu Adjuster