


# ABRASIVES

## DURA-GREEN® STONES














- Impregnated with silicon carbide
- Long lasting, fast cutting & always concentric in motion
- Use **Dura-Green Stones for:** composites, amalgams, precious alloys, non-precious alloys & porcelains
- Use **Dura-Green Fine Stones for:** composites & porcelains
- See page 24, Unmounted Finishing Wheels section, for Dura-Green Wheels

SHAPE CATEGORY	CN1	CN2	CN3	CY2	CY3	CY4	FL1	FL2	FL3	FL4	FL5
SHAPE											
ISO#	024	037	040	025	028	045	030	025	030	045	065
HP (12 PK)	PN 0002	-	-	0006	0007	0008	0009	0010	0011	0012	0013
HP FINE (12 PK)	PN -	-	-	-	-	-	-	-	011F	-	-
HP (72 PK)	PN -	-	-	-	-	-	-	HG10	-	-	-
CA (12 PK)	PN 0052	0053	0054	-	0057	-	0059	0060	0061	0062	-
FG (12 PK)	PN 0101	-	-	0103	-	-	-	0104	-	-	-

SHAPE CATEGORY	IC1	IC2	IC3	IC4	IC5	IC6	IC7	IC8	IC9	KN1	KN3
SHAPE											
ISO#	030	032	050	045	070	065	125	104	135	050	070
HP (12 PK)	PN 0014	0015	0016	0017	0018	0019	0020	0021	0022	0023	0025
HP FINE (12 PK)	PN -	-	-	-	-	-	-	-	022F	-	-
HP (72 PK)	PN HG14	-	HG16	-	HG18	-	HG20	-	HG22	-	-
CA (12 PK)	PN 0064	-	0066	-	0068	-	-	-	-	-	-
FG (12 PK)	PN -	0105	-	-	-	-	-	-	-	-	-

SHAPE CATEGORY	KN4	KN5	KN6	KN7	PC1	PC2	RD1	RD2	RD3	RD4	RE1
SHAPE											
ISO#	080	060	100	125	030	050	030	024	036	050	062
HP (12 PK)	PN 0026	0027	0028	0029	0030	0031	0032	0033	0034	0035	0036
HP FINE (12 PK)	PN -	-	-	-	030F	031F	-	-	-	-	-
HP (72 PK)	PN -	-	-	HG29	-	-	-	-	-	-	-
CA (12 PK)	PN -	-	-	-	0080	0081	0082	0083	0084	0085	0086
FG (12 PK)	PN -	-	-	-	-	-	0107	0108	-	-	-



SHAPE CATEGORY	RE2	TC1	TC2	TC3	TC4	WH1	WH2	WH3	WH4	WH5	WH6
SHAPE											
ISO#	062	020	025	030	040	040	050	065	060	065	130
HP (12 PK)	PN 0037	0039	0040	0041	0042	0043	0044	0045	0046	0047	0048
HP FINE (12 PK)	PN -	039F	-	-	042F	-	-	045F	-	-	048F
HP (72 PK)	PN -	-	-	-	HG42	-	-	HG45	-	-	HG48
CA (12 PK)	PN -	0089	0090	-	0092	-	-	0095	0096	-	0098
FG (12 PK)	PN -	0109	-	-	-	0110	-	-	-	-	-

### DURA-GREEN® STONES

**PN**

**Product**

See charts . . . . .HP Shank (12 per pack)

See charts . . . . .HP Shank (1/2 Gross/One Box)

See charts . . . . .CA & FG Shank (12 per pack)