


# ABRASIVES

## DURA-GREEN® STONES














- Impregnated with silicon carbide
- Long lasting, fast cutting & always concentric in motion
- Use **Dura-Green Stones for:** composites, amalgams, precious alloys, non-precious alloys & porcelains
- Use **Dura-Green Fine Stones for:** composites & porcelains
- See page 24, Unmounted Finishing Wheels section, for Dura-Green Wheels

| SHAPE CATEGORY  | CN1     | CN2  | CN3  | CY2  | CY3  | CY4  | FL1  | FL2  | FL3  | FL4  | FL5  |
|-----------------|---------|------|------|------|------|------|------|------|------|------|------|
| SHAPE           |         |      |      |      |      |      |      |      |      |      |      |
| ISO#            | 024     | 037  | 040  | 025  | 028  | 045  | 030  | 025  | 030  | 045  | 065  |
| HP (12 PK)      | PN 0002 | -    | -    | 0006 | 0007 | 0008 | 0009 | 0010 | 0011 | 0012 | 0013 |
| HP FINE (12 PK) | PN -    | -    | -    | -    | -    | -    | -    | -    | 011F | -    | -    |
| HP (72 PK)      | PN -    | -    | -    | -    | -    | -    | -    | HG10 | -    | -    | -    |
| CA (12 PK)      | PN 0052 | 0053 | 0054 | -    | 0057 | -    | 0059 | 0060 | 0061 | 0062 | -    |
| FG (12 PK)      | PN 0101 | -    | -    | 0103 | -    | -    | -    | 0104 | -    | -    | -    |

| SHAPE CATEGORY  | IC1     | IC2  | IC3  | IC4  | IC5  | IC6  | IC7  | IC8  | IC9  | KN1  | KN3  |
|-----------------|---------|------|------|------|------|------|------|------|------|------|------|
| SHAPE           |         |      |      |      |      |      |      |      |      |      |      |
| ISO#            | 030     | 032  | 050  | 045  | 070  | 065  | 125  | 104  | 135  | 050  | 070  |
| HP (12 PK)      | PN 0014 | 0015 | 0016 | 0017 | 0018 | 0019 | 0020 | 0021 | 0022 | 0023 | 0025 |
| HP FINE (12 PK) | PN -    | -    | -    | -    | -    | -    | -    | -    | 022F | -    | -    |
| HP (72 PK)      | PN HG14 | -    | HG16 | -    | HG18 | -    | HG20 | -    | HG22 | -    | -    |
| CA (12 PK)      | PN 0064 | -    | 0066 | -    | 0068 | -    | -    | -    | -    | -    | -    |
| FG (12 PK)      | PN -    | 0105 | -    | -    | -    | -    | -    | -    | -    | -    | -    |

| SHAPE CATEGORY  | KN4     | KN5  | KN6  | KN7  | PC1  | PC2  | RD1  | RD2  | RD3  | RD4  | RE1  |
|-----------------|---------|------|------|------|------|------|------|------|------|------|------|
| SHAPE           |         |      |      |      |      |      |      |      |      |      |      |
| ISO#            | 080     | 060  | 100  | 125  | 030  | 050  | 030  | 024  | 036  | 050  | 062  |
| HP (12 PK)      | PN 0026 | 0027 | 0028 | 0029 | 0030 | 0031 | 0032 | 0033 | 0034 | 0035 | 0036 |
| HP FINE (12 PK) | PN -    | -    | -    | -    | 030F | 031F | -    | -    | -    | -    | -    |
| HP (72 PK)      | PN -    | -    | -    | HG29 | -    | -    | -    | -    | -    | -    | -    |
| CA (12 PK)      | PN -    | -    | -    | -    | 0080 | 0081 | 0082 | 0083 | 0084 | 0085 | 0086 |
| FG (12 PK)      | PN -    | -    | -    | -    | -    | -    | 0107 | 0108 | -    | -    | -    |



| SHAPE CATEGORY  | RE2   | TC1   | TC2   | TC3   | TC4   | WH1   | WH2   | WH3   | WH4   | WH5   | WH6   |
|-----------------|---|---|---|---|---|---|---|---|---|---|---|
| SHAPE           |  |  |  |  |  |  |  |  |  |  |  |
| ISO#            | 062   | 020   | 025   | 030   | 040   | 040   | 050   | 065   | 060   | 065   | 130   |
| HP (12 PK)      | PN 0037   | 0039  | 0040  | 0041  | 0042  | 0043  | 0044  | 0045  | 0046  | 0047  | 0048  |
| HP FINE (12 PK) | PN -  | 039F  | -   | -   | 042F  | -   | -   | 045F  | -   | -   | 048F  |
| HP (72 PK)      | PN -  | -   | -   | -   | HG42  | -   | -   | HG45  | -   | -   | HG48  |
| CA (12 PK)      | PN -  | 0089  | 0090  | -   | 0092  | -   | -   | 0095  | 0096  | -   | 0098  |
| FG (12 PK)      | PN -  | 0109  | -   | -   | -   | 0110  | -   | -   | -   | -   | -   |

### DURA-GREEN® STONES

**PN**

**Product**

See charts . . . . .HP Shank (12 per pack)

See charts . . . . .HP Shank (1/2 Gross/One Box)

See charts . . . . .CA & FG Shank (12 per pack)