

GLASIONOMER BASE CEMENT (DENTIN REPLACEMENT)



- Recommended as a supportive, protective base under stress-bearing amalgam, composite, porcelain & indirect composite restorations
- Biocompatible & chemically bonds to dentin and enamel
- High fluoride release reduces the occurrence of secondary caries
- Compressive strength (147 MPa)
- High level of radiopacity
- Low solubility
- Easy mixing & handling
- Pediatric temporary cusp replacement

GLASIONOMER BASE CEMENT

PN	Product
1110	Glaslonomer Base Cement Kit 15g yellow color powder, 5ml liquid, mixing pad, spatula, dispenser & instructions
1111	Yellow Powder 15g
1113	White Powder 15g
1114	Base Cement Liquid 5ml