

HY-BOND® ZINC PHOSPHATE CEMENT



- Seals dentinal tubules & minimizes post-operative sensitivity
- Excellent for seating permanent prostheses
- Film thickness (25µm)
- High compressive strength (98 MPa)
- High level of radiopacity
- Apply directly to prepared dentin
- Neutralizes the harshness of phosphoric acid

HY-BOND® ZINC PHOSPHATE CEMENT

PN

Product

1170 Zinc Phosphate Cement Kit
60g powder, 25ml liquid, dispenser & instructions

1171 Zinc Phosphate Powder 60g

1172 Zinc Phosphate Liquid 25ml