



HY-BOND® ZINC PHOSPHATE CEMENT



- Seals dentinal tubules & minimizes post-operative sensitivity
- Excellent for seating permanent prostheses
- Film thickness (25µm)
- High compressive strength (98 MPa)
- High level of radiopacity
- Apply directly to prepared dentin
- Neutralizes the harshness of phosphoric acid

HY-BOND® ZINC PHOSPHATE CEMENT

PN

Product

- 1170** **Zinc Phosphate Cement Kit**
60g powder, 25ml liquid, dispenser & instructions
- 1171** **Zinc Phosphate Powder 60g**
- 1172** **Zinc Phosphate Liquid 25ml**