













# DIAMONDS/CARBIDES

## HYBRID POINTS® (DIAMONDS)

• Hybrid Points cut like a 6 fluted carbide bur

SHAPE CATEGORY	2	4	6	34	35	245	331	556/56	557/57	558/58
SHAPE										
ISO#	001/010	001/014	001/018	010/008	010/010	238/008	237/010	108/009	107/010	107/012
LENGTH	1.0	1.0	1.6	1.0	0.8	3.0	1.4	2.9	4.0	4.5
GRIT	R	R	R	R	R	R	R	R	R	R
PN (1 PK)	0901-1	0903-1	0905-1	0907-1	0909-1	0911-1	0913-1	0915-1	0917-1	0919-1

SHAPE CATEGORY	700/170	701/701
SHAPE		
ISO#	170/010	170/012
LENGTH	4.0	4.0
GRIT	R	R
PN (1 PK)	0921-1	0923-1

## SPECIALTY DIAMONDS

PN

Product

See charts . . . . . Hybrid Points (1 per pack)

0900 . . . . . Hybrid Points Kit

12 Hybrid Points, 10 T&F Hybrid Points