

## MONOCeM® SELF-ADHESIVE RESIN CEMENT



- Self-etch & self-adhesive
- Bonds to all substrates with high retention values
- Moisture tolerant with two second tack cure for easy clean-up
- Eliminates sensitivity & has high fluoride release
- 100% polymerization/dual cure: three minute auto cure, 20 second light cure
- Command set, curing begins at restoration's placement
- Very low film thickness (<12µm)

## MONOCeM® SELF-ADHESIVE RESIN CEMENT

### Product

PN

### **3208A . . . . . MonoCem Translucent Kit**

2 (3.5g) Translucent, 15 applicator tips & 2 OneGloss PS Midi-Point polishers

### **3209A . . . . . MonoCem Bleach White Kit**

2 (3.5g) Bleach White, 15 applicator tips & 2 OneGloss PS Midi-Point polishers

### **3211 . . . . . MonoCem Auto-mix Tips**

20 applicator tips