

ABRASIVES

ONEGLOSS® & ONEGLOSS® PS



- OneGloss aluminum oxide polishers are disposable
- OneGloss complete kit with stainless steel CA mandrels
- OneGloss PS aluminum oxide polishers are mounted on sturdy plastic mandrels
- OneGloss PS are disposable & individually wrapped
- Available in three different shapes
- Durable enough to withstand polishing an entire arch
- **Use OneGloss & OneGloss PS for:** composites, temporary crowns, removal of orthodontic cements & removal of stains from enamel

SHAPE CATEGORY		Cup	Midi-Pt.	IC	Mini-Pt.	Assort. Kit	CA Mandrel
SHAPE						See Descriptions Below	
ISO#		060	060	100	030		
ONEGLOSS							
CA (50 pk)	PN	0181	0182	0183	0186 (12 pk)	0180	0435 (6 pk)
ONEGLOSS PS							
CA (50 pk)	PN	0176	0177	0178	-	0175	-

ONEGLOSS® & ONEGLOSS® PS

PN

Product

0180 OneGloss Set

20 of each Cup, Midi-Point & IC shapes, plus 3 CA mandrels

See chart. OneGloss Re ill (50 per box)

Cup, Midi-Point or IC

0186 OneGloss Mini-Point Re ill (12 per pack)

0435 OneGloss CA Mandrel (6 per box)

0175 OneGloss PS Assortment Kit

15 Cups, 20 Midi-Points, 15 IC

See chart. OneGloss PS (50 per box)

Cup, Midi-Point or IC

0185 OneGloss/SuperBuff Combo Set

20 each of OneGloss Cup, Midi-Point & IC shapes, 25 SuperBuffs disks, 5 CA mandrels