

# OneGloss PS

## One-Step Finishing and Polishing System



### OneGloss PS One-Step Finishing and Polishing System

*Shofu is dedicated to manufacturing the highest quality products utilizing our innovative technologies to enhance minimally invasive cosmetic dentistry.*

shofu.com  
800.827.4638  
customer-service@shofu.com



### FAST AND EASY

Save chairside time when finishing and polishing all composite restorations with **OneGloss PS**. When used on the slow speed with water spray; a firmer pressure will remove composite while softer pressure will polish to a satin finish.



Finish and polish in one step.



Lasts longer than the competitors.




Does not scratch enamel.

### PROVEN QUALITY AND SUPERIOR VALUE



When compared to Enhance®, the leading competitor, **OneGloss PS** is rated much higher; finishing and polishing more restorations before showing any signs of deterioration.\* Even better, **OneGloss PS** offers 25% more product for about the same price.†



†Information is based on Henry Schein Dental Catalog.

## OneGloss

### PROVEN QUALITY

*For over 100 years, Shofu Dental's dedication to quality and value have made us one of the world's most respected dental products manufacturers.*



Scan Here  
Instructions  
for Use

SHAPE CATEGORY	Cup	Midi-Point	IC	Assort. Kit
Shape				-
<b>OneGloss</b>				
CA (50 pk) PN	0181	0182	0183	0180
<b>OneGloss PS</b>				
CA (50 pk) PN	0176	0177	0178	0175