

Pediatric Solutions

FIT SA

Self-Adhesive Flowable
Featuring Therapeutic
Giomer Technology



Anti-Bacterial
Protection




Prevents Secondary
Decay*



Acid Neutralization



Faster, Easier &
More Profitable



Bioactive Giomer
Technology**



Conventional
Composites





Mounika Falemban, DMD



Fast, easy, and reliable PRRs

- > Revolutionary self-adhesive technology saves time and money by eliminating the need for bonding agent
- > Superior bond strength for long-term durability
- > Therapeutic Giomer Technology is clinically proven to reduce secondary decay
- > Low flow viscosity for perfect placement

Resin Tag Infiltration



FIT SA F03

Shade	Item No.
A1	Y2550
A2	Y2551

CX-PLUS

Glass Ionomer
Luting Cement

- > Chemically bonds to dentin and enamel for excellent marginal seal
- > High levels of fluoride release and excellent radiopacity
- > Capsules or Powder & Liquid delivery systems
- > Value priced compared to major competitors



Fred Margolis, DDS



CX-Plus 50 Capsules

Item No. 3277

CX-Plus Liquid/Powder Kit

Item No. 1166



BeautiSealant

No-Etch Fluoride Releasing Pit and Fissure Sealant



Anti-Bacterial Protection



Prevents Secondary Decay*




Acid Neutralization




Faster, Easier & More Profitable



Bioactive Gioner Technology**



Conventional Sealants



Dr. Satoshi Fukumoto, Tohoku University, Japan†




- > No acid etching, saving time and money!
- > Superior bond strength for long-term durability
- > Therapeutic Gioner Technology is clinically proven to reduce secondary decay
- > BPA- and HEMA-free chemistry




BeautiSealant System
Item No. 1798

Deep Penetration into Fissures†



50 x 200 μm



100 x 100 μm

OneGloss PS

Single Step Finishing and Polishing Solution



- > Fast and easy one-step polishing of composite restorations
- > Firm pressure removes composite and soft pressure polishes to a "satin finish"
- > Individually wrapped and disposable to prevent cross-contamination
- > Available in three shapes for ease of use
- > Value priced compared to major competitors



OneGloss PS Assort. Kit

Item No. 0175

OneGloss PS Cup (50)

Item No. 0176

OneGloss PS Midi-Pt. (50)

Item No. 0177

OneGloss PS IC (50)

Item No. 0178



BIOACTIVE.
GIONER TECHNOLOGY.
THE BETTER CHOICE.

Visit shofu.com
Call 800.827.4638

Shofu Dental Corporation San Marcos CA 92078



*8-Year independent clinical study recorded 100% retention rate, no secondary caries, failures, or postoperative sensitivity.

**Asahi University School of Dentistry. FIT SA Dental Advisor evaluation is available at shofu.com or by request.

SN17-0721