

Polishing Kits

Finishing and Polishing Made Easy



Shofu Polishing Kits
Finishing and Polishing
Made Easy

**SERVE
PATIENTS
INCREASE
PRODUCTION
GO HOME
ON TIME!**

Shofu is dedicated to manufacturing the highest quality products utilizing our innovative technologies to enhance minimally invasive cosmetic dentistry.

shofu.com
800.827.4638
customer-service@shofu.com

SIMPLIFIED KIT SELECTION

Shofu Dental is going back to basics with a portfolio of new Abrasive Kits to allow for stress-free selection of the ideal kit regardless of the restorative material being polished. Each kit has the perfect assortment of abrasives and polishers to guarantee fast, predictable, and aesthetic results. Simply match the Abrasive Kit to the Dental material — it's that easy!



PROVEN QUALITY

With a century of experience and a tradition of excellence, Shofu is the original polishing company. This is the reason why so many Dental Schools and Clinicians rely on Shofu's legendary abrasives and polishers countless times a day. Only trust your restorations and reputation to the industry leader.



Dr. Jack Griffin



Dr. Jack Griffin



Dr. Sushil Koirala



Dr. Howard Glazer



Dr. Muneki Hirayama



SUPERIOR VALUE

You don't have to pay more to get the best! In addition to superior performance, Shofu's abrasive and polishing kits offer value that no competitor can match. Many of our competitors have reduced the number of abrasive and polishing points in their kits. Shofu kits contain an ample supply of components to polish better.

Composite Polishing Kit CA

For finishing and polishing all types of composites to a lustrous shine.

Porcelain & Enamel Polishing Kit CA

For adjusting, finishing, and polishing porcelain restoration and enamel defects.

Amalgam & Metals Polishing Kit FG

For high speed finishing, polishing and super-polishing amalgam, precious, and non-precious alloys.

Amalgam & Metals Polishing Kit CA

For slow speed finishing, polishing and super-polishing amalgam, precious, and non-precious alloys.

Acrylic Polishing Kit HP

For chairside or laboratory adjusting and polishing acrylic dentures and appliances.

Zirconia Polishing Kit CA

Diamond impregnated polishing system for polishing today's harder materials to a lustrous shine.

shofu.com
800.827.4638
customer-service@shofu.com




- COMPOSITE POLISHING KIT CA (\$79.95)***
(Item No. 0310B)
 - Four CompoSite Midi-Pts., two of each Regular & Fine
 - Four CompoSite Cups, two of each Regular & Fine
 - Two CompoSite Knives, one of each Regular & Fine
 - One Dura-White® FL2
 - One Dura-White® RD2



- PORCELAIN & ENAMEL POLISHING KIT CA (\$79.95)***
(Item No. 0317B)
 - Five Ceramisté® Midi-Pts., two of each Standard, Ultra & one Ultra II
 - Three Ceramisté® Cups, one of each Standard, Ultra & one Ultra II
 - Two Ceramisté® Wheels, one of each Standard & Ultra
 - One Dura-White® FL2
 - One Dura-White® RD2



- AMALGAM & METALS POLISHING KIT FG (\$79.95)***
(Item No. 0309B)
 - Seven Mini-Pts., three of each Brownie® & Greenie® & one Supergreenie
 - Five Cups, two of each Brownie® & Greenie® & one Supergreenie



- AMALGAM & METALS POLISHING KIT CA (\$79.95)***
(Item No. 0308B)
 - Four Mini-Pts., two of each Brownie® & Greenie®
 - Four Cups, two of each Brownie® & Greenie®
 - Two PC2's, one of each Brownie® & Greenie®
 - Two WH6's, one of each Brownie® & Greenie®



- ACRYLIC POLISHING KIT HP (\$134.95)***
(Item No. 0321B)
 - Three AcryPoint BP1's, one of each Coarse, Medium & Fine
 - Three AcryPoint PC's, two of each Coarse, Medium & Fine
 - One Robot® Carbide HP Cutter SH 79E
 - One Robot® Carbide HP Cutter SH 251E



- ZIRCONIA POLISHING KIT CA (\$179.95)***
(Item No. 0333B)
 - Three ZiLMaster Bullets, Coarse, Medium & Fine
 - Three ZiLMaster Cups, Coarse, Medium & Fine
 - Three ZiLMaster Mini-Pts., Coarse, Medium & Fine

*Suggested Retail
SN9-0522