

CEMENTS

RESICEM® UNIVERSAL RESIN CEMENT



- Designed specially for zirconia & Shofu HC Blocks – highest reported bond strengths of any resin cement
- Formulated adhesion & aesthetics
- Customize to be universal for all materials
- Self-conditioning
- High fluoride release
- Very low film thickness (9µm)

RESICEM® UNIVERSAL RESIN CEMENT

PN

Product

3220 ResiCem Basic Kit

2 (9.08g) syringes, 3ml Primer A, 2.5ml Primer B, 20 mixing tips, 1 fine microtip handle, 50 fine microtip applicators & petit dish

3221 ResiCem

9.08g syringe & 10 mixing tips

3223 ResiCem Paste Opaque (5ml)

3227 Primer A (3ml)

3228 Primer B (2.5ml)

3232 AZ Primer (5ml)

for bonding to zirconia (zirconia primer)

1740 ML Primer (5ml)

for bonding to metal (metal primer)

3234 Porcelain Primer (5ml)

for bonding to porcelain & indirect composites