

DIAMONDS/CARBIDES

SHOFU ROBOT® CARBIDE FG FINISHERS



- High quality 10 & 12 blade tungsten carbide finishing head on a stainless steel shank
- Surface finish of all types of restorative materials, both anterior & posterior regions
- **Use Robot Carbide FG Finishers for:** all restorative materials

SHAPE CATEGORY	SH 41	SH 46	SH 46	SH 379	SH 379	SH 379	SH 375R
SHAPE							
HEAD SIZE (MM)	014	014	018	014	018	023	014
HEAD LENGTH (MM)	-	3.5	3.8	3.1	3.5	4.2	8.0
# OF BLADES	12	10	10	10	10	10	10
MAXIMUM RPM	450,000	450,000	450,000	450,000	450,000	300,000	300,000
FG (3 pk)	PN 0703	0709	0710	0717	0718	0719	0727

SHAPE CATEGORY	SH 247	SH 246	SH 246	SH 246
SHAPE				
HEAD SIZE (MM)	009	009	010	012
HEAD LENGTH (MM)	3.2	3.6	3.6	3.6
# OF BLADES	10	10	10	10
MAXIMUM RPM	450,000	450,000	450,000	450,000
FG (3 pk)	PN 0736	0741	0742	0743

PRICING - SHOFU ROBOT® CARBIDE FG FINISHERS

PN

Product

See chartsFG Finishers (3 per pack)