

DIAMONDS/CARBIDES

SHOFU ROBOT® CARBIDE FG FINISHERS



- High quality 10 & 12 blade tungsten carbide finishing head on a stainless steel shank
- Surface finish of all types of restorative materials, both anterior & posterior regions
- **Use Robot Carbide FG Finishers for:** all restorative materials

SHAPE CATEGORY	SH 41	SH 41	SH 46	SH 46	SH 379	SH 379	SH 379
SHAPE							
HEAD SIZE (MM)	014	018	014	018	014	018	023
HEAD LENGTH (MM)	-	-	3.5	3.8	3.1	3.5	4.2
# OF BLADES	12	12	10	10	10	10	10
MAXIMUM RPM	450,000	450,000	450,000	450,000	450,000	450,000	300,000
FG (3 PK)	PN 0703	0704	0709	0710	0717	0718	0719

SHAPE CATEGORY	SH 375R	SH 247	SH 246	SH 246	SH 246
SHAPE					
HEAD SIZE (MM)	014	009	009	010	012
HEAD LENGTH (MM)	8.0	3.2	3.6	3.6	3.6
# OF BLADES	10	10	10	10	10
MAXIMUM RPM	300,000	450,000	450,000	450,000	450,000
FG (3 PK)	PN 0727	0736	0741	0742	0743

SHOFU ROBOT® CARBIDE FG FINISHERS

PN Product
See chartsFG Finishers (3 per pack)