

DIAMONDS/CARBIDES

SHOFU ROBOT® CARBIDE HP CUTTERS



- High quality tungsten carbide finishing head on a stainless steel shank
- Universal cross cut
- Use **Robot Carbide HP Cutters for:** predominantly base, noble & high noble metal alloys, chrome cobalt, model casts, plaster & acrylic

SHAPE CATEGORY	SH 71E	SH 71E	SH 72E	SH 73E	SH 73E	SH 77E	SH 77E	SH 78E	SH 79E	SH 79E
SHAPE										
HEAD SIZE (MM)	023	050	060	014	060	023	060	023	040	050
HEAD LENGTH (MM)	-	-	12.0	3.1	10.0	5.0	11.0	6.0	13.0	13.0
# OF BLADES	10/9	16/10	18/12	10/6	18/12	14/12	18/12	14/9	12/8	16/13
MAXIMUM RPM	50,000	50,000	50,000	50,000	50,000	50,000	50,000	50,000	50,000	50,000
HP (1/PK)	PN 0602	0603	0605	0606	0609	0610	0613	0614	0616	0617

SHAPE CATEGORY	SH 251E	SH 351E	SH 137E	SH 261E	SH 295E	SH 129E	SH 33
SHAPE							
HEAD SIZE (MM)	060	060	016	023	023	023	009
HEAD LENGTH (MM)	14.0	11.0	4.0	13.0	15.0	8.0	4.2
# OF BLADES	18/12	18/12	10/6	10/9	14/12	14/12	6
MAXIMUM RPM	50,000	50,000	50,000	50,000	50,000	50,000	50,000
HP (1/PK)	PN 0620	0623	0626	0630	0631	0634	0667(6PK)

SHOFU ROBOT® CARBIDE HP CUTTERS

PN

Product

See chartsHP Cutters (1 per pack)

0602, 0606, 0610, 0614, 0626, 0630, 0631, 0634

See chartsHP Cutters (1 per pack)

0603, 0605, 0609, 0613, 0605, 0616, 0617, 0620, 0623

0667HP Cutters (6 per pack)