

Solidilite V & Sublite V

PURE AESTHETICS



LIGHT CURING UNITS



Visit www.shofu.com or call 800.827.4638

HIGH PERFORMANCE LIGHT TECHNOLOGY OPTIMIZED FOR CERAMAGE® AND SOLIDEX®

Maximize the strength and working efficiency of your lab processed indirect composite restorations and light cured acrylics with the Solidilite V light curing unit and the easy to use Sublite V tacking unit.



SOLIDILITE V

The New Shofu Solidilite V features a wide polymerization chamber and easy to use presets. The presets are for one minute P1, three minutes P2 and five minutes P3. This light curing unit has four strategically placed powerful JCR 110v 150W/S halogen lights and a rotating turntable for fast, efficient, and shadow free polymerization. Solidilite V can accommodate various voltage situations (100 to 240v, 50-60Hz) with the turn of a simple voltage selection button in the back of the unit.



SUBLITE V

The Sublite V is designed for initial flash polymerization with easy one-touch operation. The one touch hand operated switch irradiates for five seconds and by continuously pressing the switch the irradiation will last for 20 seconds, then automatically turn off. The compact size (W4½xD5xH10 inches) allows the unit to be conveniently placed on the bench.

FEATURES AND BENEFITS

Compatible voltage	Accommodates 100-240v (50-60 Hz)
Multiple presets (1 min., 3 min. and 5 min.)	One button start
Manual override setting ability	Establish your own presets
Shadow-free illumination	Fast and complete polymerization
7- segment LED display	Visible display of time remaining
Manual up/down control of curing table	Easy and accurate adjustments
Lamp life indicator	Detects when lamps are dead
Easy lamp replacement	Saves time and reduces down time
Solidilite's footprint (W9xD9¾xH13½ in.)	Easily fits in work space
Adjustable Sublite light guard	Safely positions for easy working
Sublite's footprint (W4½xD5xH10 in.)	Easily fits in work space



SOLIDILITE V (PN 5105)
SUBLITE V (PN 5107)



SHOFU INC. 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983, Japan (www.shofu.co.jp)
SHOFU INC. Singapore Branch 10 Science Park Road, #03-12 The Alpha, Science Park II, Singapore 117684 (www.shofu.com.sg)
SHOFU DENTAL TRADING (SHANGHAI) CO., LTD. No. 645 Jiye Road, Sheshan Industrial Park, Songjiang, Shanghai 201602, China (www.shofu.com.cn)
SHOFU DENTAL CORPORATION 1225 Stone Drive, San Marcos, CA 92078-4059, USA (www.shofu.com)
SHOFU UK Riverside House, River Lawn Road, Tonbridge, Kent, TN9 1EP, UK (www.shofu.co.uk)
SHOFU DENTAL GmbH [CE] [REP] Am Brüll 17, 40878 Ratingen, Germany (www.shofu.de)
 SHOFU is a registered trademark of SHOFU INC. All other trademarks and registered trademarks are the property of their respective holders. SHOFU INC. reserves the right to change specifications without notice.
 © Copyright 2001 SHOFU INC. All rights reserved.

SHOFU



CERAMICS