

Super-Snap

Finishing and Polishing Disk System



Super-Snap Finishing and Polishing Disk System

Shofu is dedicated to manufacturing the highest quality products utilizing our innovative technologies to enhance minimally invasive cosmetic dentistry.

shofu.com
800.827.4638
customer-service@shofu.com



SUPERIOR DESIGN

Most major competitors feature a metal ring attachment system, which may damage the surface of a newly placed restoration. **Super-Snap's** superior design offers a completely flat polishing surface that will never harm a restoration, while also providing more surface area for faster polishing.



Super-Snap disks feature superior design



Super-Snap Disks are flexible and ultra-thin

PROVEN QUALITY AND SUPERIOR VALUE



Super-Snap sets the industry standard for finishing and polishing anterior composite restorations. Available in two diameters, simply go from the coarse grit to the fine grit in order to achieve a glossy finish. Even better, Super-Snap is 30% less expensive than the Sof-Lex™ disc polishing system.*



*Information is based on Henry Schein Dental Catalog.

Super-Snap



Scan Here
Instructions
for Use

- ☐ **Super-Snap Rainbow Kit**
(Item No. 0500)
 - 180 assorted disks (100 standard & 80 mini)
 - 40 Polystrips
 - 2 Dura-White Stones (CN1 & FL2)
 - 1 Composite Fine Midi-Point
 - 4 CA Metal Mandrels



- ☐ **Super-Snap Singles Standard Size Disk (12mm)**
(Item No. 0502)
 - 1 Disk of Coarse, Medium, Fine & Super Fine grits
 - Super-Snap plastic CA mandrel
- ☐ **Super-Snap Singles Mini Size Disk (8mm)**
(Item No. 0503)
 - 1 Disk of Coarse, Medium, Fine & Super Fine grits
 - Super-Snap plastic CA mandrel

SHAPE CATEGORY	Contouring		Finishing				Polishing		Super-Polishing	
Super-Snap Shape										
Color/Grit	Black (coarse)		Violet (medium)				Green (fine)		Red (super fine)	
Refills (50 pk) PN	L506	L507	L508	L528	L509	L522	L501	L519	L502	L521