



SUPER-SNAP® - FINISHING & POLISHING



- Elastic silicone shank mount for easy placing of the disk onto the mandrel
- Manufactured without a metal center
- Very flexible ultra-thin disks allow easy access to interproximal areas
- Two sizes – four grits (silicon carbide & aluminum oxide)
- **Use Super-Snap for:** composites

SUPER-SNAP® - FINISHING & POLISHING

PN

Product

0500 Super-Snap Rainbow Technique Kit

180 assorted disks (100 Standard & 80 Mini), 40 Polystrips, 2 Dura-White Stones (CN1 & FL2), 1 CompoSite Fine Midi-Point, 4 CA metal mandrels arranged in a semicircular rainbow organizer tray

0505 Super-Snap Mini Kit

48 assorted disks (32 Standard & 16 Mini), 1 CompoSite Fine Midi-Point, 1 CA metal mandrel

0501 Super-Snap Singles

Mini & Standard - Coarse, Medium, Fine & Super Fine grits, Super-Snap plastic CA mandrel (8 disks per Single - 48 Singles per box)

0502 Super-Snap Singles Standard Size Disk (12mm)

1 disk of Coarse, Medium, Fine & Super Fine grits, Super-Snap plastic CA mandrel (4 disks per Single - 48 Singles per box)











0503 Super-Snap Singles Mini Size Disk (8mm)

1 disk of Coarse, Medium, Fine & Super Fine grits, Super-Snap plastic CA mandrel (4 disks per Single - 48 Singles per box)

ABRASIVES

SUPER-SNAP® DISK REPLACEMENTS

- Black (Coarse): Two sizes – safe side down, disposable silicon carbide disk
- Violet (Medium): Two sizes – safe side down, double sided & safe side up, disposable silicon carbide disk
- Green (Fine): Two sizes – double sided, aluminum oxide grit, safe for enamel
- Red (Super Fine): Two sizes – double sided, aluminum oxide grit, safe for enamel

SHAPE CATEGORY	Contouring		Finishing				Polishing		Super-Polishing		
SHAPE		MINI 				MINI 		MINI 		MINI 	
COLOR/GRIT	Black (Coarse)		Violet (Medium)				Green (Fine)		Red (Super Fine)		
REFILLS (50 PK)	PN	L506	L507	L508	L528	L509	L522	L501	L519	L502	L521

SUPER-SNAP® DISK REPLACEMENTS

PN

Product

See chart **Standard & Mini (50 per box)**