











DIAMONDS/CARBIDES

T&F HYBRID POINTS® (DIAMONDS)

• T&F Hybrid Points cut like a 12 fluted carbide bur

SHAPE CATEGORY	7006	7104	7106	7404	7406	7803	7901	7903	LT1	LT2
SHAPE										
ISO#	001/018	254/014	254/018	277/014	277/018	499/012	243/009	243/012	164/013	165/011
LENGTH	1.6	3.0	3.5	3.0	3.5	4.0	3.7	4.2	6.0	8.0
GRIT	SF	SF	SF	SF	SF	SF	SF	SF	SF	SF
PN (1 PK)	0925-1	0927-1	0929-1	0931-1	0933-1	0935-1	0937-1	0939-1	0941-1	0943-1

T&F HYBRID POINTS®

PN

Product

See charts T&F Hybrid Points (1 per pack)

0950 T&F Hybrid Points Kit (all shapes)