

Shade	W1	W2	W3	A1	A2	A3	A3.5	A4	rootA
<b>Opaque Liner</b>	W1O	W2O	W3O	A1O	A2O	A3O	A3.5O	A4O	rootAO
<b>Cervical</b>	—	—	—	—	—	A3B:2 AC:1	A3.5B:1 AC:1	A4B:1 AC:1	AC
<b>Margin</b>	W1M	W1M:2 NM:1	W1M:1 NM:2	NM	NM:1 A3M:1	A3M	A3M:2 rootAM:1	A3M:1 rootAM:2	rootAM
<b>Opaque Dentin</b>	OD-W1	OD-W1:2 OD-N:1	OD-W1:1 OD-N:2	OD-N	OD-N:1 OD-A3:1	OD-A3	OD-A3:2 OD-rootA:1	OD-A3:1 OD-rootA:2	OD-rootA
<b>Body</b>	W1B	W2B	W3B	A1B	A2B	A3B	A3.5B	A4B	rootAB
<b>Enamel</b>	Opal 56	Opal 56:2 Opal 57:1	Opal 56:1 Opal 57:2	Opal 57	Opal 58	Opal 59	Opal 59:1 Opal 60:1	Opal 60	

Shade	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
<b>Opaque Liner</b>	B1O	B2O	B3O	B4O	C1O	C2O	C3O	C4O	D2O	D3O	D4O
<b>Cervical</b>	—	—	B3B:1 BC:1	BC	—	C2B:2 CC:1	C3B:1 CC:1	CC	DC		
<b>Margin</b>	NM:1 B2M:1	B2M	B2M:1 B4M:1	B4M	NM:1 C2M:1	C2M	C2M:1 C4M:1	C4M	NM:1 D3M:1	D3M	D3M:1 B4M:1
<b>Opaque Dentin</b>	OD-N:1 OD-B2:1	OD-B2	OD-B2:1 OD-B4:1	OD-B4	OD-N:1 OD-C2:1	OD-C2	OD-C2:1 OD-C4:1	OD-C4	OD-N:1 OD-D3:1	OD-D3	OD-D3:1 OD-B4:1
<b>Body</b>	B1B	B2B	B3B	B4B	C1B	C2B	C3B	C4B	D2B	D3B	D4B
<b>Enamel</b>	Opal 57	Opal 58	Opal 59	Opal 60	Opal 58	Opal 59	Opal 60	Opal 58	Opal 59		

<b>Opal Effect</b>	Opal T	Opal SL	Opal WE	Opal MI	Opal OC	Opal AM-R	Opal AM-Y	Opal AM-V
--------------------	--------	---------	---------	---------	---------	-----------	-----------	-----------

<b>Enamel Effect</b>	BT	OT	PT	GT	T-GLASS	<b>Cervical Trans</b>	CT-CL	CT-W	CT-A	CT-B	CT-R
----------------------	----	----	----	----	---------	-----------------------	-------	------	------	------	------

<b>Color Effect</b>	MP	MY	Mlv	RED	Y	O	G	W	<b>Gum Color</b>	Gum-LP	Gum-DP
	OM-Y	OM-LP	OM-DP	LPM	CLM	<b>Correction</b>			ADD-ON B	ADD-ON T	CPM-FINE

※印刷物のため実物とは色調が異なる場合があります。