










ZiLMASTER POINTS



- ZiLMaster is a three-step system for contouring, finishing, pre-polishing & polishing
- Diamond impregnated silicone polishers that come in three different grits
- Coarse is suitable for abrasive prepolishing & finishing of areas or zones that still need minor material reduction
- Medium is recommended for contouring & smoothing anatomical details & preparing surfaces for high-gloss polishing
- Fine is designed for efficient high-gloss polishing without the additional use of polishing pastes
- **Use ZiLMaster for:** zirconia & lithium disilicate

SHAPE CATEGORY	Bullet	Knife	Cup	Mini-Pt.	KN7	WH6	Unm. Wheel	Kits
SHAPE								See Assort. Below
ISO#	040	090	065	030	125	120	170	
COARSE (green band)								
HP (3 PK)	PN 0640	-	-	0641	0642	0643	0645 (1 per pack)	0644
CA (3 PK)	PN 0658	0662	0663	0661	-	-		0664
MEDIUM (blue band)								
HP (3 PK)	PN 0646	-	-	0647	0648	0649	0651 (1 per pack)	0650
CA (3 PK)	PN 0665	0659	0668	0666	-	-		0669
FINE (yellow band)								
HP (3 PK)	PN 0652	-	-	0653	0654	0655	0657 (1 per pack)	0656
CA (3 PK)	PN 0670	0672	0673	0671	-	-		0674

ZiLMASTER POINTS

PN

Product

0644 **Coarse Assortment Package - HP Shank**

Bullet, Mini-Point, KN7, WH6

0650 **Medium Assortment Package - HP Shank**

Bullet, Mini-Point, KN7, WH6

0656 **Fine Assortment Package - HP Shank**

Bullet, Mini-Point, KN7, WH6

ABRASIVES

ZiLMASTER POINTS

ZiLMASTER POINTS

PN	Product
0664	Coarse Assortment Package - CA Shank Bullet, Knife, Cup, Mini-Point
0669	Medium Assortment Package - CA Shank Bullet, Knife, Cup, Mini-Point
0674	Fine Assortment Package - CA Shank Bullet, Knife, Cup, Mini-Point
See chart	CA & HP Shank (3 per pack)
See chart	Unmounted Wheel (1 per pack)